

# Detectable Changes in the Level of Cerebral Activation are Related to Auditory Working Memory Load of the Frontal Lobe: A Functional Near Infrared Spectroscopy Study.

G.T. Voelbel<sup>1,2</sup>, J. Lengenfelder<sup>1,2</sup>, G. Wylie<sup>1,2</sup>, R. L. Barbour<sup>3,4</sup>, Y. Pei<sup>3</sup>, A. Smith<sup>1</sup>, & J. DeLuca<sup>1,2</sup>

<sup>1</sup>Kessler Medical Rehabilitation Research & Education Center, West Orange, NJ <sup>2</sup>University of Medicine and Dentistry of New Jersey

<sup>3</sup>NIRx Medical Technologies <sup>4</sup>SUNY Downstate Medical Center

## Introduction

- > This study examined the concentration of oxygenated hemoglobin (oxyHb) during the N-Back task, a verbal working memory task, with a functional near infrared spectroscopy system (See Figure 1).
- > The N-Back paradigm used here has 4 levels of working memory cognitive load ("0-back," "1-back," "2-back," and "3-back," see Figure 2).
- > The relationship between the level of working memory load and hemodynamic response of oxyHb was explored.

Figure 1. Multi-channel Continuous Wave Near Infrared Imager and Helmet




Figure 2. N-Back Paradigm Conditions





Figure 3. Three dimensional model of source-detector pairs of the frontal lobe.



## Methods

### Participants

- > 9 Right-Handed Healthy Adults
- > Age: 30.7 (11.3) years
- > Education: 15.9 (2.7) years
- > Free of substance abuse and major psychiatric disorders
- > No history of neurological disease or trauma

### Behavioral Task

#### N-Back Task

- > Consonant letters were presented every three seconds auditorily in 24 second epochs.
- > Four Conditions (0-back, 1-back, 2-back, 3-back)
- > Each condition randomly presented 3 times

### Procedures

#### Apparatus

- > Multi-channel continuous wave near infrared imager (NIRx Medical Technologies; see Figure 2)
- > 30 source and 30 detector optodes (900 channels)
- > Dual wavelength of near infrared light (760nm and 830nm)
- > Optodes placed on forehead 10% above nasion in a 10 cm by 3 cm rectangle configuration, Figure 3.

### Data Preprocessing and Analysis

#### Preprocessing

- > Near Infrared Analysis, Visualization and Imaging (NAVI) software (NIRx Medical Technologies, LLC)
- > Low-band pass filter (.15 Hz)
- > 15% mean Coefficients of Variation threshold
- > Oxy-Hb concentration was modeled with a modified Lambert-Beer analysis for each time point in each voxel of modeled space (Figure 3).
- > Data converted to Analyze format

#### Data Analysis

- > AFNI image analysis software
- > Time-series deconvolved for each N-back Condition
- > N-Back conditions compared across participants with t-tests (random effects analysis).
- > Results corrected for multiple comparisons
  - >  $\alpha = 0.05$ , cluster size = 31 contiguous voxels.
- > Plots of oxyHb
  - > Data baseline corrected over first 7 sec
  - > Then averaged across subjects
  - > Plots shown at locations of reliable difference

Funded by the National Institute of Neurological Disorders and Stroke (1F32 NS055509-01 & R41 NS050007) and the Henry H Kessler Foundation. Presented at the Cognitive Neuroscience Society Annual Meeting, 2007

## Results

Figure 4. 1-Back Minus 0-Back Tasks

Area of activation: BA 47, Right Middle Frontal Gyrus




Figure 5. 2-Back Minus 1-Back Conditions

Area of activation: BA 10/46, Right Middle Frontal Gyrus  
BA 47, Left Middle Frontal Gyrus.





Figure 6. 3-Back Minus 2-Back Conditions

Area of Activation: BA 46, Right Inferior Frontal Gyrus



## Conclusions

- > Increased oxy-Hb is associated with greater verbal working memory cognitive demand primarily in right ventro-lateral prefrontal cortex (middle/inferior frontal gyrus).